

Aquinity² P35/70

Produces pure water & ultra pure water




**WATER
PURIFIER**






The **Aquinity² P35/70** system is designed for the production of pure water ($< 1 \mu\text{S}/\text{cm}$) and ultra pure water ($0.055 \mu\text{S}/\text{cm}$) out of tap water.

The background of the slide features a dynamic splash of clear water against a white background. The water is captured in mid-air, with numerous droplets and a large, frothy plume on the right side. The splash creates a sense of freshness and purity, which is thematically appropriate for a presentation about water filtration technology.

The Aquinity² P35/70 system is equipped with a reverse osmosis (RO) and conditioning module to produce deionised water (DI water) into the external storage reservoir (35 L or 70 L) with a production rate of 10 L/h and an access to tap DI water directly. Out of the storage tank the ultra pure water is produced with mixed bed cartridges and then can be tapped with a flexible dispenser. During non-use automatic recirculation within the Aquinity² system ensures a permanent high quality of ultrapure water. All components of the Aquinity² are assembled in an especially designed housing. The system can be easily opened by removing the side panel to access the ultra pure module and UV lamp (185 nm/ 254 nm).



	Tap Water
	< 10 µS/ cm
	< 1 µS/ cm, Type II
	0.055 µS/ cm, Type I

Technical Specifications




conductivity	< 1 $\mu\text{S}/\text{cm}$; ASTM Type II + 0.055 $\mu\text{S}/\text{cm}$; ASTM Type I
resistivity	> 1 $\text{M}\Omega\cdot\text{cm}$; ASTM Type II + 18.2 $\text{M}\Omega\cdot\text{cm}$; ASTM Type I
total organic carbon (TOC)	< 10 ppb Reagent < 5 ppb Life Science < 3 ppb Analytical
dispensing flow rate	2.0 L/min (1.5 L/min Life Science)
productivity rate	10 L/h, optionally 20 L/h
bacteria	< 0.01 cfu/mL*
particulate	> 0.2 μm , less than 1 particulate/mL
pyrogen (endotoxins)	< 0.001 EU/mL
RNAse DNAse	< 0.5 pg/mL* < 10 pg/mL*
dimensions, weight, power	650 x 560 ⁺ x 540 ⁺ (W x D x H), 16-20 kg, 110-230 V

* with Life Science model or microbiological endfilter

⁺ with dispenser 650 x 590 x 750 mm

Feed Water Requirements

type of feed water	tap water
feed water conductivity	< 1400 $\mu\text{S}/\text{cm}$
inlet pressure	1.5 to 6 bar
free chlorine	< 0.1 mg/L
Silt Density Index (SDI)	< 3
pH	3 to 9
temperature	5 to 25 °C

A large, dynamic splash of water is visible in the background, extending across the bottom half of the page. The water is captured in mid-air, creating a sense of movement and freshness. The splash is centered and spreads outwards, with many small droplets visible. The overall color palette is light blue and white, matching the clean, technical theme of the document.

Configurations

Aquinity² P35/70 is available in different configurations to fit the specific requirements of pure and ultrapure water quality for different applications.



Aquinity ² (10 L/h)	UV reactor	UF module	TOC monitoring	Cat. No. (35 L reservoir)	Cat. No. (70 L reservoir)	Aquinity ² 20 (20 L/h)
Reagent	-	-	-	114-0050	114-0060	-
Analytical	+	-	-	114-0051	114-0061	-
Life Science	+	+	-	114-0052	114-0062	-
Analytical TI	+	-	+	114-0056	114-0066	114-0074
Life Science TI	+	+	+	114-0057	114-0067	114-0075

Options



TOC monitoring (TI version)

The TOC monitoring during production and intermittent measurements during non-use periods allows to check the organic content in water continuously between 1 and 999 ppb.




volumetric dispensing

Our dispenser allows the volumetric controlled dispensing of water with an increment of 0.1 L and a tap volume from 0.1 to 99 L. The system prevents overflow of containers and allows to dispense water without supervising.


wall mounted configuration & bench integrated (BI version)

To save valuable benchtop space, the system is available in two main space-saving configurations. Either the entire system can be wall-mounted and an inclined display for a more convenient operation or ordered as a bench- integrated configuration. Alternatively, only the system or the tank can be fitted in this way. This results in a wide variety of combinations. Furthermore it is possible to install the display and dispenser on the wall.




User Interface & Software


Currently recorded data and warning messages will be displayed on the touch screen monitor. The software furthermore supports the user with a self-diagnostic module which reduced service time and costs.




Helping tools to maintain system



Trackable values up to 1 year



Warning of customer in case of over limit and service case



Main Screen: The Software allows the user to see all information to use maintenance tool and track back historical values up to 1 year.





Consumables & Spare Parts

Description	Cat. No.
Conditioning Module P	190-0086
Ultra pure Module MemPak Life Science/Analytical	190-0087/ 190-0088
vent filter for external tank	190-0085
UV lamp (185 nm/ 254 nm)	921-0138
Submersible UV lamp for reservoir	921-0659
Pretreatment Module ProPak R10	290-0065
Ultrafiltration Module	190-0096
Final Filter, capsule 0.2 µm	190-0013*
Disinfection Tablets	290-0227

* 190-0305 Final Filter to reduce Endotoxin, DNA + RNase

Made in Germany

Your trusted partner for more than 30 years



membraPure GmbH
Wolfgang-Küntscher-Str. 14
16761 Hennigsdorf, Germany
Tel.: +49 3302 - 201 20 0



info@membrapure.de
www.membrapure.de